Instructions for use

(EN)	Instructions for use	3
(DE)	Einbauanleitung	5
(FR)	Guide de fabrication	7
(ES)	Instrucciones para el uso	9
(IT)	Istruzioni per l'uso	11
(DA)	Brugsanvisning	13
(SV)	Bruksanvisning	15
(NL)	Gebruiksaanwijzing voor de vervaadiging	17
(PO)	Instruções de Utilização	19
(JA)	取扱説明書	21

Junior 4-Prong Socket Adapter

Product Number : AJ-114040

Life Without Limitations

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For secure fastening use the loctite recommended and apply rated torque twice.

Instructions for integration into laminated socket.

Note: For more complex positions, the socket can be laminated in two layers. The pronged adapter should be shaped and fixed to the inner layer with putty. All surfaces should be roughened prior to the second lamination.

When shaping the pronged adapter, care must be taken to prevent disortion of the distal threads.

- I. Isolate cast with PVA bag
- II. Apply Double layer of Perlon Stockinet
- III. Apply single 50mm carbon tapes crossed over distal end of cast, in AP/ML direction
- IV. Shape Prong Adapter and position on distal end of

- model in correct alignment and rotation
- V. Screw Lamination cap into distal threads
 - VI. Apply clay to Pinch Bolt and Pinch Bolt housing to protect from resin
 - VII. Tie 50mm carbon tapes into grooves on Adapter, position in AP/ML positions
 - VIII. Lay up material and combinations of material should be appropriate to the individuals
 - IX. Layers of material may be trimmed to different lengths to provide a thinner lighter proximal socket
 - X. Tie 50mm carbon tapes into grooves on Adapter, position in AP/ML positions
 - XI. Apply double layer of Perlon Stockinet
 - XII. Laminate with 80:20 acrylic lamination resin

The Lamination Cap should be removed following lamination by unscrewing from Prong Adapter.

- 2 year limited warranty

Steel to ISO10328, load level A60, rated to 100lbs / 45kg, for normal use.

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